





260-BURKIT01 HP Bur FG



All burs are 44.5 mm and are used in Low-speed 1:1 handpieces and straight nose cone. Suited for larger and exotic animals. To avoid overheating water is advisable during use.




Round burs – To be used for caries excavation, cavity preparation, root canal access and osseous surgery.


Burs	#	Quantity	Description
	260-H111012	2	Round, size 3, 44.5 mm
	260-H111023	2	Round, size 8, 44.5 mm

Cross-cut taper fissure burs – Ideal for sectioning multi-rooted teeth and reducing crown height.

Burs	#	Quantity	Description
	260-H3111557	2	Flat-End Cross-Cut Fissure, # 557, 44.5 mm
	260-H3311021	2	Flat-End Cross-Cut Fissure, # 703L, 44.5 mm

Inverted cone – With a slightly rounded edge on blade corners, which produces a smoother and superior cut. Generally used for entering carious pits, producing an undercut in cavity preparations with a flat floor.

Burs	#	Quantity	Description
	260-H211018	1	Inverted Cone, # 39, 44.5 mm

Burs	#	Quantity	Description
	260-H79G11050	1	Coarse Standard Cutter with Spiral Blade Tothing, 44.5 mm



D I S P O M E D